

Category	Cat. No.	Product	Description	Type	Size	Enzyme	Application
Human	<b>Virus</b>						
	IP10010	iNNOPLEX™ Enterovirus VP1 Detection Kit(1st)	Enterovirus VP1 gene	RT-PCR(1st)	96T	Maxime RT-PCR	Human
	IP10020	iNNOPLEX™ Enterovirus VP1 Detection Kit(2nd)	Enterovirus VP1 gene	nested-PCR(2nd)	96T	Maxime i-StarTaq	Human
	IP10030	iNNOPLEX™ Enterovirus NCR Detection Kit	Enterovirus 5'NCR gene	RT-PCR	96T	Maxime RT-PCR	Human
	IP10040	iNNOPLEX™ Rota VP4/VP7 Detection Kit	Rotavirus VP4, VP7 gene	RT-PCR, duplex	96T	Maxime RT-PCR	Human
	D80100	iNNOPLEX™ West Nile Detection Kit(1st : RT-PCR)	West Nile virus	RT-PCR	96T	Maxime RT-PCR	Human
	D80200	iNNOPLEX™ West Nile Detection Kit(2nd : PCR)	West Nile virus	PCR	96T	Maxime i-StarTaq	Human
	D80110	iNNOPLEX™ JEV Detection Kit(1st : RT-PCR)	Japanese Encephalitis virus	RT-PCR	96T	Maxime RT-PCR	Human
	D80120	iNNOPLEX™ JEV Detection Kit(2nd : PCR)	Japanese Encephalitis virus	PCR	96T	Maxime i-D Taq	Human
	<b>Parasite</b>						
	IP10090	iNNOPLEX™ TSUTSUJ Detection Kit	<i>Tsutsugamushi</i>	nested-PCR	48T	i-D Maxime	Human
	<b>Food pathogen</b>						
	D10370	VeTeK™ LEPTO Detection Kit	<i>Leptospira spp.</i>	PCR	96T	Maxime i-StarTaq	Human
	D10380	VeTeK™ SAL Detection Kit	<i>Salmonella spp.</i>	PCR	96T	Maxime i-StarTaq	Human
D10390	VeTeK™ CAMPY Detection Kit	<i>Campylobacter jejuni</i>	PCR	96T	Maxime i-StarTaq	Human	
D10400	VeTeK™ CLOS Detection Kit	<i>Clostridium perfringens</i>	PCR	96T	Maxime i-StarTaq	Human	
D50050	VeTeK™ VIBRIO Detection Kit	<i>Vibrio spp.</i>	PCR	96T	Maxime i-StarTaq	Human	
	<b>Virus</b>						
	D10230	VeTeK™ CPV Detection Kit	CPV type-1 and type-2	PCR, duplex	96T	Maxime i-StarTaq	Canine
	D10240	VeTeK™ CAdV Detection Kit	Canine adeno virus	PCR	96T	Maxime i-StarTaq	Canine
	D10250	VeTeK™ CDV Detection Kit	Canine distemper virus	Nested-PCR	48T	Maxime RT-PCR + i-StarTaq	Canine
	D10260	VeTeK™ CCoV Detection Kit	Canine corona virus	RT-PCR	96T	Maxime RT-PCR	Canine
	D10270	VeTeK™ CPIV Detection Kit	Canine parainfluenza virus	RT-PCR	96T	Maxime RT-PCR	Canine
	D10310	VeTeK™ RV Detection Kit	Rabies virus	RT-PCR	96T	Maxime i-StarTaq	Zoonosis
	D10330	VeTeK™ CHV Detection Kit	Canine Herpesvirus	PCR	96T	Maxime i-StarTaq	Canine
	D20110	VeTeK™ FIV Detection Kit	Feline immunodeficiency virus	PCR	96T	Maxime i-StarTaq	Feline
	D20120	VeTeK™ FeLV Detection Kit	Feline leukemia virus	RT-PCR	96T	Maxime RT-PCR	Feline
	D20140	VeTeK™ FPV Detection Kit	Feline panleukopenia	PCR	96T	Maxime i-StarTaq	Feline
	D20150	VeTeK™ FCoV Detection Kit	FCoV or FIPV	RT-PCR	96T	Maxime i-StarTaq	Feline
	D20180	VeTeK™ FHV Detection Kit	Feline herpes virus	PCR	96T	Maxime i-StarTaq	Feline
	D20190	VeTeK™ FCV Detection Kit	Feline calicivirus	RT-PCR	96T	Maxime i-StarTaq	Feline
	D30160	VeTeK™ CAV Detection Kit	Chicken anemia virus	PCR	96T	Maxime i-StarTaq	Avian
	D30180	VeTeK™ ALV Detection Kit	Avian leukosis virus	PCR	96T	Maxime i-StarTaq	Avian
	D30200	VeTeK™ MDV Detection Kit	Marek disease virus	Nested-PCR	48T	Maxime i-StarTaq	Avian
	D30210	VeTeK™ ILTV Detection Kit	Infectious laryngotracheitis virus	PCR	96T	Maxime i-StarTaq	Avian
	D30220	VeTeK™ NDV Detection Kit	NDV(common & patho)	RT-PCR	48T	Maxime i-StarTaq	Avian
	D30230	VeTeK™ IBV Detection Kit	Infectious bronchitis virus	RT-PCR	96T	Maxime i-StarTaq	Avian
	D30240	VeTeK™ AIV Detection Kit		RT-PCR(duplex)	96T	Maxime RT-PCR	Avian
	D30250	VeTeK™ AIV H5 Detection Kit	AIV serotype H5	RT-PCR	96T	Maxime RT-PCR	Avian
	D30260	VeTeK™ AIV H7 Detection Kit	AIV serotype H7	RT-PCR	96T	Maxime RT-PCR	Avian
	D30270	VeTeK™ AIV H9 Detection Kit	AIV serotype H9	RT-PCR	96T	Maxime RT-PCR	Avian
	D30290	VeTeK™ IBDV Detection Kit	Infectious bursal disease virus	RT-PCR	96T	Maxime RT-PCR	Avian
	D40170	VeTeK™ PCV Detection Kit	Procine circo virus	PCR	96T	Maxime i-StarTaq	Porcine
	D40180	VeTeK™ PRRS Detection Kit	PRRSV	RT-PCR	96T	Maxime i-StarTaq	Porcine
	D40190	VeTeK™ SIV Detection Kit	Swine influenza virus	RT-PCR	96T	Maxime RT-PCR	Porcine
	D40200	VeTeK™ TGE/PED Detection Kit	TGEV, PEDV	RT-PCR, duplex	96T	Maxime i-StarTaq	Porcine
	D40220	VeTeK™ PCV-2 Detection Kit	Porcine Circovirus type 2	PCR	96T	Maxime i-StarTaq	Porcine
	D40230	VeTeK™ ADV Detection Kit	Aujeszky's disease virus	PCR	96T	Maxime i-StarTaq	Porcine
	D60040	VeTeK™ BRV Detection Kit	Bovine Rotavirus	RT-PCR	96T	Maxime i-StarTaq	Bovine
	D60050	VeTeK™ BCoV Detection Kit	Bovine Coronavirus	RT-PCR	96T	Maxime i-StarTaq	Bovine
	D60060	VeTeK™ BVDV Detection Kit	BVDV type-1,2, BDV, CSFV	RT-PCR	96T	Maxime i-StarTaq	Ruminant, Porcine
	D60070	VeTeK™ BD-Multi Detection Kit	BRV, BCoV, BVDV type-1,2, BDV, CSFV	RT-PCR, multiplex	50T	Maxime i-StarTaq	Ruminant, Porcine

Category	Cat. No.	Product	Description	Type	Size	Enzyme	Application
Animal	D60130	iNNOPLEX™ IBR Detection Kit	Bovine Herpesvirus type I	PCR	96T	Maxime i-StarTaq	Bovine
	D80100	iNNOPLEX™ West Nile Detection Kit(1st : RT-PCR)	West Nile virus	RT-PCR	96T	Maxime RT-PCR	Animal
	D80200	iNNOPLEX™ West Nile Detection Kit(2nd : PCR)	West Nile virus	PCR	96T	Maxime i-StarTaq	Animal
	D80110	iNNOPLEX™ JEV Detection Kit(1st : RT-PCR)	Japanese Encephalitis virus	RT-PCR	96T	Maxime RT-PCR	Animal
	D80120	iNNOPLEX™ JEV Detection Kit(2nd : PCR)	Japanese Encephalitis virus	PCR	96T	Maxime i-D Taq	Animal
	D90020	iNNOPLEX™ Akabane Detection Kit(1st : RT-PCR)	Akabane virus	RT-PCR(1st)	96T	Maxime RT-PCR	Ruminant
	D90120	iNNOPLEX™ Akabane Detection Kit(2nd : PCR)	Akabane virus	nested-PCR(2nd)	96T	Maxime i-StarTaq	Ruminant
	D90130	iNNOPLEX™ Aino Detection Kit(1st : RT-PCR)	Aino virus	RT-PCR(1st)	96T	Maxime RT-PCR	Ruminant
	D90050	iNNOPLEX™ Aino Detection Kit(2nd : PCR)	Aino virus	nested-PCR(2nd)	96T	Maxime i-StarTaq	Ruminant
	D90060	iNNOPLEX™ Chuzan Detection Kit(1st : RT-PCR)	Chuzan virus	RT-PCR(1st)	96T	Maxime RT-PCR	Ruminant
	D90070	iNNOPLEX™ Chuzan Detection Kit(2nd : PCR)	Chuzan virus	nested-PCR(2nd)	96T	Maxime i-StarTaq	Ruminant
	D90080	iNNOPLEX™ BEFV Detection Kit(1st : RT-PCR)	Bovine ephemeral fever virus	RT-PCR(1st)	96T	Maxime RT-PCR	Ruminant
	D90090	iNNOPLEX™ BEFV Detection Kit(2nd : PCR)	Bovine ephemeral fever virus	nested-PCR(2nd)	96T	Maxime i-StarTaq	Ruminant
	D90100	iNNOPLEX™ Ibaraki Detection Kit(1st : RT-PCR)	Ibaraki virus	RT-PCR(1st)	96T	Maxime RT-PCR	Ruminant
	D90110	iNNOPLEX™ Ibaraki Detection Kit	Ibaraki virus	nested-PCR(2nd)	96T	Maxime i-StarTaq	Ruminant
	D90140	iNNOPLEX™ Simbu Detection Kit	Simbu serogroup	RT-PCR	96T	Maxime RT-PCR	Ruminant
	D90150	iNNOPLEX™ EHDV Detection Kit	EHDV(Palyam) serogroup	RT-PCR	96T	Maxime RT-PCR	Ruminant
	D90180	iNNOPLEX™ RHD Detection Kit	Rabbit haemorrhagic disease virus	RT-PCR	96T	Maxime RT-PCR	Rabbit
<b>Bacteria</b>							
	D10280	VeTeK™ BAB-C Detection Kit	<i>Babesia canis &amp; B. gibsoni</i>	PCR	96T	Maxime i-StarTaq	Canine
	D10290	VeTeK™ HBN-C Detection Kit	<i>Mycoplasma haemocanis</i>	PCR	96T	Maxime i-StarTaq	Canine
	D10300	VeTeK™ EHR Detection Kit	<i>Ehrlichia spp.</i>	PCR	96T	Maxime i-StarTaq	Zoonosis
	D10350	VeTeK™ Myco Detection Kit	<i>Mycoplasma spp.</i>	PCR	96T	Maxime i-StarTaq	Animal
	D10370	VeTeK™ LEPTO Detection Kit	<i>Leptospira spp.</i>	PCR	96T	Maxime i-StarTaq	Animal
	D10380	VeTeK™ SAL Detection Kit	<i>Salmonella spp.</i>	PCR	96T	Maxime i-StarTaq	Animal
	D10390	VeTeK™ CAMPY Detection Kit	<i>Campylobacter jejuni</i>	PCR	96T	Maxime i-StarTaq	Animal
	D10400	VeTeK™ CLOS Detection Kit	<i>Clostridium perfringens</i>	PCR	96T	Maxime i-StarTaq	Animal
	D10410	VeTeK™ RICK Detection Kit	<i>Rickettsia rickettsii</i>	PCR	96T	Maxime i-StarTaq	Zoonosis
	D10420	VeTeK™ BRU Detection Kit	<i>Brucella spp.</i>	PCR	96T	Maxime i-StarTaq	Zoonosis
	D20130	VeTeK™ HBN-F Detection Kit	<i>Mycoplasma haemofelis</i>	PCR	96T	Maxime i-StarTaq	Feline
	D20170	VeTeK™ CHLA Detection Kit	<i>Chlamydia psittaci</i>	PCR	96T	Maxime i-StarTaq	Feline, Avian
	D20200	VeTeK™ TOXO Detection Kit	<i>Toxoplasma gondii</i>	PCR	96T	Maxime i-StarTaq	Zoonosis
	D30170	VeTeK™ APEC Detection Kit	Avian pathogenic <i>Escherichia coli</i>	PCR	96T	Maxime i-StarTaq	Avian
	D30330	VeTeK™ MS/MG Detection Kit	<i>M. gallisepticum &amp; M. synoviae</i>	Nested-PCR	24T	Maxime i-StarTaq	Avian
	D40130	VeTeK™ PAS-A Detection Kit	<i>Pasteurella multocida</i> type A	PCR	96T	Maxime i-StarTaq	Animal
	D40140	VeTeK™ PAS-D Detection Kit	<i>Pasteurella multocida</i> type D	PCR	96T	Maxime i-StarTaq	Animal
	D40150	VeTeK™ BOR Detection Kit	<i>Bordetella bronchiseptica</i>	PCR	96T	Maxime i-StarTaq	Animal
	D40160	VeTeK™ Myco-P Detection Kit	<i>Mycoplasma hyopneumoniae</i>	PCR	96T	Maxime i-StarTaq	Animal
	D40240	VeTeK™ PAS Detection Kit	<i>Pasteurella multocida</i>	PCR	96T	Maxime i-StarTaq	Animal
	D50050	VeTeK™ VIBRIO Detection Kit	<i>Vibrio spp.</i>	PCR	96T	Maxime i-StarTaq	Animal
	D60080	iNNOPLEX™ Blackleg Detection Kit	<i>Clostridium chauvoei</i>	PCR	96T	Maxime i-StarTaq	Animal
	D60090	iNNOPLEX™ Anthrax Detection Kit	<i>Bacillus anthracis</i>	PCR(duplex)	96T	Maxime i-StarTaq	Animal, Bio-terror
	D60100	iNNOPLEX™ TB Detection Kit(IS6110, IS1081)	<i>M. Tuberculosis</i> (IS6110, IS1081)	PCR(duplex)	96T	Maxime i-StarTaq	Animal
	D60110	iNNOPLEX™ TB Detection Kit(16S)	<i>M. Tuberculosis</i> (16S)	PCR	96T	Maxime i-StarTaq	Animal
<b>Parasite</b>							
	D10320	VeTeK™ HW Detection Kit	<i>Diro laria immitis</i>	PCR	96T	Maxime i-StarTaq	Canin, Feline
	D10340	VeTeK™ GIA Detection Kit	<i>Giardia spp.</i>	PCR	96T	Maxime i-StarTaq	Animal
	D60140	iNNOPLEX™ Neospora Detection Kit	<i>Neospora caninum</i>	PCR	96T	Maxime i-StarTaq	Animal
<b>Fungi</b>							
	D10360	VeTeK™ ASP Detection Kit	<i>Aspergillus spp.</i>	PCR	96T	Maxime i-StarTaq	Animal
Fish	<b>Virus</b>						
	D50060	VeTeK™ WSSV Detection kit (ONE-STEP Type)	White spot syndrome virus	PCR	96T	Maxime i-StarTaq	Shrimp
	D50061	VeTeK™ WSSV Detection kit (Nested PCR Type)	White spot syndrome virus	nested-PCR	48T	Maxime i-StarTaq	Shrimp

Category	Cat. No.	Product	Description	Type	Size	Enzyme	Application	
Fish	D50070	iNNOPLEX™ IHNV Detection Kit	IHNV	PCR			Shrimp	
	<b>Bacteria</b>							
	D50050	VeTeK™ VIBRIO Detection Kit	Vibrio spp.	PCR	96T	Maxime i-StarTaq	Fish, shrimp	
Plant	<b>Virus</b>							
	DP10010	i-PSTVd Detection Kit	Potato Spindle Tuber Viroid	RT-PCR	96 T	i-D Maxime RT-PCR	Potato	
	DP10011	i-PVS Detection Kit	Potato Virus S	RT-PCR	96 T	Maxime RT-PCR	Potato	
	DP10012	i-PVY Detection Kit	Potato Virus Y	RT-PCR	96 T	Maxime RT-PCR	Potato	
	DP10013	i-PLRV Detection Kit	Potato Leaf Roll Virus	RT-PCR	96 T	Maxime RT-PCR	Potato	
	DP10014	i-PVX Detection Kit	Potato Virus X	RT-PCR	96 T	Maxime RT-PCR	Potato	
	DP10015	i-PVM Detection Kit	Potato Virus M	RT-PCR	96 T	Maxime RT-PCR	Potato	
	DP10021	i-RSV™ Detection Kit	Rice Strip Tenuivirus	RT-PCR	96 T	Maxime RT-PCR	Rice	
	<b>Bacteria</b>							
		DP10016	iNNOPLEX™ PBW R1 Detection Kit	Potato bacterial wilt	PCR	96 T	Maxime i-StarTaq	Potato
		DP10017	iNNOPLEX™ PBW R3 Detection Kit	Potato bacterial wilt	PCR	96 T	Maxime i-StarTaq	Potato
	<b>Prep</b>							
		17411	Plant Pathogenic Viral dsRNA Isolation Kit	Double strand RNA isolation(resin type)	RNA Prep	50 rxn		Potato
	17412	Double-RNA Viral dsRNA Extraction Kit(for Plant Tissue)	Double strand RNA isolation(membrane type)	RNA Prep	50 rxn		Plant	
	DP10020	i-RSV™ Isolation Kit	Rice Strip Tenuivirus	RNA Prep	96 rxn		Rice	